

Brecknell

HS-300

Portable Physician Scale



Simplicity of design and rugged construction are the key features of the HS-300 Portable Physician Scale. This scale provides accurate weight in lb/kg on easy-to-read LCD backlit display (batteries, adapter and COMM cable included).

It has a large weighing range and can measure height from 2' to 7' x 0.125" / 60 cm to 210 cm x 0.1 cm.

Constructed of mild steel with a durable powder coated white finish, this scale is suitable for almost any application requiring accuracy and medium to heavy weight capacity.



Accuracy

Within 0.1% of full load



Capacity

660 lb x 0.2 lb / 5 oz / 300 kg x 0.1 kg



Power

9 VDC, 600mA power adapter (included)



Display

17 mm (0.625") high backlit LCD



Construction

Mild steel with a durable powder coated white finish



Operator Keys

Hold/Print, Units, Tare, On/Off/Zero

VALUE IN
WEIGHING

HS-300 Portable Physician Scale

Specifications

Readability: Accurate full range weights (lb/kg) on backlit LCD display. Built in Weight Hold feature. Height gauge to 120 cm / 7'

Features: Rollers included for easy transportation. Powder coated steel construction with durable finish. Easy to operate digital indicator with RS-232 port. Large 20.5" x 20.5" weighing platform

Platform Dimension (L x W):
20.5" x 20.5" / 521 mm x 521 mm

Handrail Height (W x H):
16" x 49" / 406 mm x 1245 mm

Scale Capacity & Height Rod Range:
660 lb x 0.2 lb / 300 kg x 0.1 kg
2' to 7' x 0.125" / 60 cm to 210 cm x 0.1 cm

Display: 17 mm (0.625") high backlit LCD

Indicator Dimension (L x W x H):
6.5" x 1.2 x 3.2" / 165 mm x 30 mm x 81 mm

Serial Port: One full duplex RS-232 with cable

Function Keys: Hold/Print, Units, Tare, On/Off/Zero

Unit: lb, kg, lb/oz

Power: 9 VDC, 600mA power adapter or four AAA batteries (both included)

Battery Life: Approximately 20 hours continuous

Operating Temperature:
41° F to 95° F / 5° C to 35° C

RS-232: Supplied complete with 10' cable (9 pin to 9 pin type connector, bi-directional full duplex)

Package Dimension (L x W x H):
30.5" x 54.5" x 9" / 775 mm x 1384 mm x 230 mm

Shipping Weight: 54 lb / 24 kg

Certifications: Full FCC compliance, ETL power adapter, Prop 65, CE, UKCA, WEEE, RoHS

Packaging Contents: Platform, Indicator, Handrail, Manual, Height Gauge, RS-232 Cable, Power Supply and Miscellaneous Hardware



Brecknell®

brecknellscales.com

brecknellscales.co.uk

Part of Avery Weigh-Tronix

Brecknell is part of Avery Weigh-Tronix. Avery Weigh-Tronix is a trademark of Illinois Tool Works group of companies whose ultimate parent company is Illinois Tool Works Inc ("Illinois Tool Works"). Copyright © 2022 Illinois Tool Works. All rights reserved. This publication is issued to provide outline information only and may not be regarded as a representation relating to the products or services concerned. This publication was correct at the time of going to print, however Avery Weigh-Tronix reserves the right to alter without notice the specification, design, price or conditions of supply of any product or service at any time.

HS-300_L_500789.indd
V2 AWT35-500789