

Brecknell

MS140-300

Portable Handrail Scale



Simple, durable, and designed with comfort in mind - the MS140-300 portable physician scale combines reliable performance with user-friendly features.

The large 20.5" x 20.5" low-profile platform and sturdy handrails with padded grips provide added stability and support, making it ideal for patients requiring extra assistance during weighing.

The steel-housed indicator features a large, easy-to-read LCD display with backlight, along with an RS232 port for efficient data transfer. A convenient Hold function freezes the last weight reading on the screen, allowing for hands-free operation.

Integrated casters make the scale easy to move when needed, while dual power options — battery or AC adapter — ensure maximum flexibility for use in a variety of healthcare settings.



Accuracy

Within 0.1% of full load



Capacity

660 lb x 0.2 lb, 300 kg x 0.1 kg,
660 lb x 5 oz



Power

Four AAA batteries or 9 VDC, 600mA
wall mount adapter (both included)



Display

Five digit 0.625" / 17 mm digits
with seven segment LCD



Construction

Stainless steel weighing surface with
mild steel foam covered handrails.
Integrated wheel kit for easy
portability.



Operator Keys

On/Off/Zero, Tare, Unit,
Hold/Print

VALUE IN
WEIGHING

MS140-300 Portable Electronic Handrail Scale

Specifications

Unit: lb, lb/oz, kg

Serial Port: One full duplex RS-232 with cable

RS-232: Supplied complete with 10' cable (9 pin to 9 pin type connector, bi-directional full duplex)

Integrated Wheel Kit: For easy portability

Weight Hold: May be used to freeze and store the displayed weight so that you may attend to the patient's needs

Operating Temperature: 41° F to 95° F / 5° C to 35° C

Scale Dimension (L x W x H):

20.5" x 20.5" x 44.25"

520 mm x 520 mm x 1135 mm

Platform Dimension (L x W x H):

20.5" x 20.5" x 2.5"

520 mm x 520 mm x 65 mm

Indicator Dimension (L x W x H):

6.5" x 1.2 x 3.2"

165 mm x 30 mm x 81 mm

Product Weight: 42 lb / 19 kg

Packaging Dimension (L x W x H):

46" x 25" x 6"

1168 mm x 635 mm x 152 mm

Shipping Weight: 55 lb / 25 kg

Certifications: FCC compliance, ETL power adapter, Prop 65

Packaging Contents: Manual, platform, three piece handles, power supply, RS-232 cable, indicator, three hex wrenches and pull through string

Warranty: 3-year limited



STAY
CONNECTED
WITH US!



@brecknell-scales



@brecknell_scales

Brecknell®

brecknellscales.com

brecknellscales.co.uk

Part of Avery Weigh-Tronix

Brecknell is part of Avery Weigh-Tronix. Avery Weigh-Tronix is a trademark of Illinois Tool Works group of companies whose ultimate parent company is Illinois Tool Works Inc ("Illinois Tool Works"). Copyright © 2025 Illinois Tool Works. All rights reserved. This publication is issued to provide outline information only and may not be regarded as a representation relating to the products or services concerned. This publication was correct at the time of going to print, however Avery Weigh-Tronix reserves the right to alter without notice the specification, design, price or conditions of supply of any product or service at any time.

MS140-300_L_500792.indd
V5 AWT35-500792