

Brecknell

MS20 Series

Medical/Veterinary Scale



The MS20 Series Medical/Veterinary Scale is designed to precisely weigh infants or small animals - ideal for applications where bench space and ergonomics matter.

The easy-to-read display connects via an extendable 16"-26" / 406 mm-660 mm coiled cable and includes a convenient Hold function that freezes the last weight shown on the screen, allowing for hands-free recording.

Operates on a single 9V battery or AC adapter, suitable for both indicator wall mounted or bench-top use.

The MS20 comes with a plastic weighing tray and the MS20S comes with a durable stainless steel weighing tray for easy cleaning and long-lasting performance.



Accuracy

Within 0.1% of full load



Capacity

44 lb x 1/2 oz / 20 kg x 10 g



Power

9V 100mA AC adapter or 9V alkaline battery (not included)



Display

Bright, large and visible 1" / 25mm LCD. The MS20 plastic model also features an audio buzzer.



Construction

Stainless steel or plastic weighing tray. Plastic housed indicator with 16"-26" / 406 mm-660 mm coiled cable.



Operator Keys

On/Off/Tare, Units, Hold

VALUE IN
WEIGHING

MS20 Series Medical/Veterinary Scale

Specifications

Capacity: 44 lb x ½ oz / 20 kg x 10 g

Display: 1" / 25 mm LCD display with four active digits

Operator Keys: On/Off/Tare, Unit, Hold

Power: 9V 100mA AC adapter or 9V alkaline battery
(not included)

Tray: Stainless steel or plastic

Units: lb, lb/oz, kg

Operating Temperature:
41° F to 95° F / 5° C to 35° C

Tray Dimension (L x W x H):
21" x 10" x 2"
533 mm x 254 mm x 51 mm

Inside Tray Dimension (L x W x H):

MS20: 20½" x 10" x 2⅝"
521 mm x 254 mm x 60 mm

MS20S: 20½" x 11½" x 2¼"
521 mm x 292 mm x 57 mm

Packaging Dimension (L x W x H):

5⅝" x 12⅜" x 10⅞"
143 mm x 313 mm x 275 mm

Shipping Weight:

MS20: 9.1 lb / 4.14 kg

MS20S: 10.75 lb / 4.89 kg

Certifications: Full FCC compliance; UL/CUL
power adapter

Warranty: 3-year limited



STAY
CONNECTED
WITH US!



@brecknell-scales



@brecknell_scales

Brecknell®

brecknellscales.com

Part of Avery Weigh-Tronix

Brecknell is part of Avery Weigh-Tronix. Avery Weigh-Tronix is a trademark of Illinois Tool Works group of companies whose ultimate parent company is Illinois Tool Works Inc ("Illinois Tool Works"). Copyright © 2025 Illinois Tool Works. All rights reserved. This publication is issued to provide outline information only and may not be regarded as a representation relating to the products or services concerned. This publication was correct at the time of going to print, however Avery Weigh-Tronix reserves the right to alter without notice the specification, design, price or conditions of supply of any product or service at any time.

MS20Series_L_09849-0014.indd
V5 09849-0014