


Brecknell

HS-200M

Portable Physician Scale



Simplicity of design and rugged construction are the key features of the HS-200M Portable Physician Scale. The scale provides accurate weight or height measurements on easy to read linear beams in bold marking. It has a large weighing range and can measure height from 2' to 7' x 0.125" / 60 cm to 210 cm x 0.1 cm. Constructed of steel with a durable white finish, this scale will provide years of service.



Construction

Steel with durable white finish



Capacity

440 lb x 4 oz / 200 kg x 100 g

VALUE IN
WEIGHING

HS-200M Portable Physician Scale

Specifications

Capacity: 440 lb x 4 oz / 200 kg x 100 g

Readability: Accurate full range weights that are dual marked in lb/kg die-cast beam. Height gauge to 7' / 210 cm

Features: Rollers included for easy transportation.
Powder coated steel construction with durable finish.
Large 14.75" x 10.8" weighing platform.

Height Rod Range:
2' to 7' x 0.125" / 60 cm to 210 cm x 0.1 cm

Platform Size (L x W):

14.75" x 10.80" / 375 mm x 275 mm

Overall Dimension (L x W x H):

22" x 10.8" x 58.5" / 559 mm x 275 mm x 1485 mm

Product Weight: 38 lb / 17 kg

Package Dimension (L x W x H):

59" x 20" x 6" / 1499 mm x 508 mm x 152 mm

Shipping Weight: 40 lb / 18 kg

Packaging Contents: Manual, platform, height scale, weighing tower, wheels and misc hardware



Brecknell®

brecknellscales.com

brecknellscales.co.uk

— *Part of Avery Weigh-Tronix*

Brecknell is part of Avery Weigh-Tronix. Avery Weigh-Tronix is a trademark of Illinois Tool Works group of companies whose ultimate parent company is Illinois Tool Works Inc ("Illinois Tool Works"). Copyright © 2023 Illinois Tool Works. All rights reserved. This publication is issued to provide outline information only and may not be regarded as a representation relating to the products or services concerned. This publication was correct at the time of going to print, however Avery Weigh-Tronix reserves the right to alter without notice the specification, design, price or conditions of supply of any product or service at any time.

HS-200M_L_500789.indd
V3 AWT35-500789