

Brecknell

MS20 Series

Medical/Veterinary Scale



The MS20 scale is particularly suited to the medical or veterinary industry.

With a capacity of 44 lb / 20 kg and readability of 1/2 oz / 10 g, the MS20 will accurately weigh infants or small animals.

The scale has a "HOLD" feature that will capture an unstable weight and keep it on the display for recording after the weight is taken off.

The MS20 is offered in two models and can operate on battery (not included) or AC adapter.

The MS20 comes with a plastic weighing tray and the MS20S comes with a stainless steel weighing tray.



Accuracy

Within 0.1% of full load



Capacity

44 lb x 1/2 oz / 20 kg x 10 g



Power

9V 100mA AC adapter or 9V alkaline battery (not included)



Display

Bright, large and visible 1" LCD



Construction

Stainless steel or plastic weighing tray.
Plastic housed indicator with 16"-26" coiled cable.



Operator Keys

On/Off/Tare, Units, Hold

VALUE IN
WEIGHING

MS20 Series Medical/Veterinary Scale

Specifications

Capacity: 44 lb x ½ oz / 20 kg x 10 g

Display: 1.0" LCD display with four active digits

Keyboard: On/Off/Tare, Unit, Hold

Power: 9V 100mA AC adapter or 9V alkaline battery
(not included)

Tray: Stainless steel or plastic

Units: lb, lb/oz, kg

Operating Temperature:
41° F to 95° F / 5° C to 35° C

Tray Dimension (L x W x H):
21" x 10" x 2"
533 mm x 254 mm x 51 mm

Inside Tray Dimension (L x W x H):

MS20: 20½" x 10" x 2⅝"
521 mm x 254 mm x 60 mm

MS20S: 20½" x 11½" x 2¼"
521 mm x 292 mm x 57 mm

Packaging Dimension (L x W x H):

5⅝" x 12⅜" x 10⅞"
143 mm x 313 mm x 275 mm

Shipping Weight:

MS20: 9.1 lb / 4.14 kg

MS20S: 10.75 lb / 4.89 kg

Certifications: Full FCC compliance; UL/CUL
power adapter

Warranty: 3-year limited



Brecknell®

brecknellscales.com

brecknellscales.co.uk



Brecknell is part of Avery Weigh-Tronix. Avery Weigh-Tronix is a trademark of Illinois Tool Works group of companies whose ultimate parent company is Illinois Tool Works Inc ("Illinois Tool Works"). Copyright © 2023 Illinois Tool Works. All rights reserved. This publication is issued to provide outline information only and may not be regarded as a representation relating to the products or services concerned. This publication was correct at the time of going to print, however Avery Weigh-Tronix reserves the right to alter without notice the specification, design, price or conditions of supply of any product or service at any time.

Part of Avery Weigh-Tronix

MS20Series_L_09849-0014.indd
V3 09849-0014