

Brecknell

MS140-300

Portable Handrail Scale



Simplicity of design and rugged construction are the key features of the MS140-300 Portable Medical Scale. Featuring a 20.5" x 20.5" / 520 mm x 520 mm, low-profile 2.5" H / 65 mm (H) platform with sturdy handrails and padded grips.

This medical scale has the ability to print weight data offering a real time saving device for doctor's offices, long-term care facilities, and other medical institutions.

Heavy-duty, yet lightweight & portable this scale is suitable for almost any application requiring accuracy and medium to heavy weight capacity. The repeatability is unbeatable and will offer years of continued usage. It comes with an easy-to-read LCD backlit display and operates on AC power. This medical scale is durable and reliable at an affordable price.



Accuracy

Within 0.1% of full load



Capacity

660 lb x 0.2 lb, 300 kg x 0.1 kg,
660 lb x 5 oz



Power

Four AAA batteries or 9 VDC, 600mA
wall mount adapter (both included)



Display

Five digit 0.625" / 17 mm digits
with seven segment LCD



Construction

Stainless steel weighing surface with
mild steel foam covered handrails.
Integrated wheel kit for easy
portability.



Operator Keys

On/Off/Zero, Tare, Unit,
Hold/Print

VALUE IN
WEIGHING

MS140-300 Portable Electronic Handrail Scale

Specifications

Unit: lb, lb/oz, kg

Serial Port: One full duplex RS-232 with cable

RS-232: Supplied complete with 10' cable (9 pin to 9 pin type connector, bi-directional full duplex)

Integrated Wheel Kit: For easy portability

Weight Hold: May be used to freeze and store the displayed weight so that you may attend to the patient's needs

Operating Temperature: 41° F to 95° F / 5° C to 35° C

Scale Dimension (L x W x H):

20.5" x 20.5" x 44.25"

520 mm x 520 mm x 1135 mm

Platform Dimension (L x W x H):

20.5" x 20.5" x 2.5"

520 mm x 520 mm x 65 mm

Indicator Dimension (L x W x H):

6.5" x 1.2 x 3.2"

165 mm x 30 mm x 81 mm

Product Weight: 42 lb / 19 kg

Packaging Dimension (L x W x H):

46" x 25" x 6"

1168 mm x 635 mm x 152 mm

Shipping Weight: 55 lb / 25 kg

Certifications: FCC compliance, ETL power adapter, Prop 65

Packaging Contents: Manual, platform, three piece handles, power supply, RS-232 cable, indicator, three hex wrenches and pull through string

Warranty: 3-year limited



Brecknell®

brecknellscales.com

brecknellscales.co.uk

Part of Avery Weigh-Tronix

Brecknell is part of Avery Weigh-Tronix. Avery Weigh-Tronix is a trademark of Illinois Tool Works group of companies whose ultimate parent company is Illinois Tool Works Inc ("Illinois Tool Works:"). Copyright © 2023 Illinois Tool Works. All rights reserved. This publication is issued to provide outline information only and may not be regarded as a representation relating to the products or services concerned. This publication was correct at the time of going to print, however Avery Weigh-Tronix reserves the right to alter without notice the specification, design, price or conditions of supply of any product or service at any time.

MS140-300_L_500792.indd
V4 AWT35-500792