




## A Range of Veterinary Scales for Small and Large Animals







Cost-effective professional veterinary scales for general purpose weighing of pets and small animals.

## Small Animals



**WS15:** The scale has a “HOLD” feature that will capture an unstable weight and

keep it on the display for recording after the weight is taken off.

- **Capacity:** 15 kg x 5 g / 30 lb x 0.01 lb
- **Construction:** Low profile stainless steel platform with mild steel understructure, stainless steel/ABS indicator will wall bracket and desk mount included
- **Platform Dimension:** 380 mm x 300 mm x 27 mm
- **Power:** Four AAA batteries or 9 VDC, 600mA adapter (both included)



**WS120:** The scale has a “HOLD” feature that will capture an unstable weight

and keep it on the display for recording after the weight is taken off.

**WS60:** The scale has a “HOLD” feature that will capture an unstable

weight and keep it on the display for recording after the weight is taken off.



- **Capacity:** 60 kg x 20 g / 130 lb x 0.05 lb
- **Construction:** Low profile stainless steel platform with mild steel understructure, stainless steel/ABS indicator will wall bracket and desk mount included
- **Platform Dimension:** 380 mm x 300 mm x 27 mm
- **Power:** Four AAA batteries or 9 VDC, 600mA adapter (both included)

- **Capacity:** 120 kg x 50 g / 260 lb x 0.1 lb
- **Construction:** Low profile stainless steel platform with mild steel understructure, stainless steel/ABS indicator will wall bracket and desk mount included.
- **Inside Tray Dimension:** 380 mm x 300 mm x 27 mm
- **Power:** Four AAA batteries or 9 VDC, 600mA adapter (both included)



Ideal for accurately weighing large animals, these scales help to remove the guesswork where treatments or administration of feed needs to be closely monitored.

## Large Animals



**WS300:** This well designed low profile scale incorporates four full load bridge load cells with self leveling, anti-slip feet to give a solid base in use.

- **Capacity:** 300 kg x 50 g
- **Construction:** Low profile stainless steel platform with mild steel understructure, stainless steel/ABS indicator will wall bracket and desk mount included
- **Platform Dimension:**  
165 mm x 30 mm x 81 mm / 6.5" x 1.2" x 3.2"
- **Power:** Four AAA batteries or 9 VDC, 600mA adapter (both included)



**PS2000:** This scale is accurate, durable, reliable, portable and cost effective putting you at ease knowing you're in good hands.

- **Capacity:** 1000 kg x 0.5 kg / 2000 lb x 1.0 lb
- **Construction:** Rugged steel with tread plate surface and stainless steel/ABS indicator
- **Platform Dimension:**  
1498 mm x 762 mm x 64 mm  
59" x 30" x 2.5"
- **Power:** 9 VDC, 600mA adapter (included) or four AAA batteries (not included)



**PS1000:** This scale is accurate, durable, reliable, portable and cost effective putting you at ease knowing you're in good hands.

- **Capacity:** 500 kg x 0.2 kg / 1000 lb x 0.5 lb
- **Construction:** Rugged steel with tread plate surface and stainless steel/ABS indicator
- **Platform Dimension:**  
1416 mm x 514 mm x 64 mm / 55.75" x 20.25" x 2.5"
- **Power:** 9 VDC, 600mA adapter (included) or four AAA batteries (not included)



**PS3000-LCD:** This scale comes with a rubber mat that reduces noise and risk of animals/goods slipping while weighing.

- **Capacity:** 1500 kg x 0.5 kg / 3000 lb x 1.0 lb
- **Construction:** Mild steel construction with safety tread plate surface
- **Platform Dimension:**  
2000 mm x 1000 mm x 100 mm  
78.75" x 39.37" x 3.93"
- **Power:** 9V 600mA universal power supply with USA/UK/EU adapters (included) or six AA batteries (not included)





## Laboratory Scales

**MBS Series:** With its fine feature set and simple functionality, this precision lab balance provides accurate weighing and counting solutions.

• **Capacity:**

300 g x 0.005 g, 600 g x 0.01 g,  
1200 g x 0.02 g, 3000 g x 0.05 g, 6000 g x 0.1 g

• **Construction:** Robust plastic housing, stainless steel top plate

• **Platform Dimension:**

300 to 600 g: 4.65" x 0.4" / 120 mm x 9.5 mm  
1200 to 6000 g: 6.9" x 6.9" x 0.3" / 175 mm x 175 mm x 8 mm

• **Power:**

AA batteries (included) with auto-off (not included)



## Hanging Scales

**ElectroSamson:** This accurate hand-held scale featured with "Hold" mode and memory tare function is ideal for poultry weighing, seed and animal feed preparation.

• **Capacity:** 10 kg x 0.01 kg / 22 lb x 0.02 lb/¼ oz  
25 kg x 0.02 kg / 55 lb x 0.05 lb/½ oz  
45 kg x 0.05 kg / 99 lb x 0.1 lb/1 oz

• **Construction:** Tough ABS case with built-in hand grip for easy fitting

• **Product Dimension:**  
11" x 5.5" x 3" / 279 mm x 140 mm x 76 mm

• **Power:** 9V alkaline battery (included)

