

 **Brecknell**

**CS-200M**

Chair Scale



User-friendly and supportive mechanical chair scale, designed for weighing individuals who have difficulty standing.

It features a durable, powder-coated steel frame and an oversized seat for maximum comfort. The armrests lift away for easy access, while the easy-to-read die-cast mechanical beam runs along the back of the chair for convenient operation.

Rear locking wheels ensure stability and safety during use, while front swivel wheels allow for smooth portability. The seat is made from easy-to-clean plastic, making the scale practical for everyday use in healthcare or assisted living environments.



**Construction**

Powder coated steel



**Capacity**

440 lb x 4 oz / 200 kg x 100 g



**Display**

Accurate full range weights that are dual marked in lb/oz and kg/g die-cast beam

VALUE IN  
WEIGHING

# CS-200M Chair Scale

## Specifications

**Capacity:** 440 lb x 4 oz / 200 kg x 100 g

**Readability:** Accurate full range weights dual marked in lb/kg die-cast beam

**Features:** Locking rear wheels for safety. Swivel front wheels for transportation. Lift way armrests for easy access.

**Seat Dimension (L x W x H):**

16" x 14" x 13"

405 mm x 350 mm x 340 mm

**Product Weight:** 70 lb / 32 kg

**Shipping Weight:** 80 lb / 36 kg

**Package Dimension (L x W x H):**

30" x 24" x 25"

762 mm x 609 mm x 635 mm

**Packaging Contents:** Manual, footrest, two armrests, scale with wheel assembly, seat, back handle, weighing tower and miscellaneous hardware

**Warranty:** 3-year limited



STAY  
CONNECTED  
WITH US!



@brecknell-scales



@brecknell\_scales

**Brecknell®**

[brecknellscales.com](http://brecknellscales.com)

[brecknellscales.co.uk](http://brecknellscales.co.uk)

Part of Avery Weigh-Tronix

Brecknell is part of Avery Weigh-Tronix. Avery Weigh-Tronix is a trademark of Illinois Tool Works group of companies whose ultimate parent company is Illinois Tool Works Inc ("Illinois Tool Works"). Copyright © 2025 Illinois Tool Works. All rights reserved. This publication is issued to provide outline information only and may not be regarded as a representation relating to the products or services concerned. This publication was correct at the time of going to print, however Avery Weigh-Tronix reserves the right to alter without notice the specification, design, price or conditions of supply of any product or service at any time.

CS-200M\_L\_500787.indd  
V4 AWT35-500787