




VALUE IN WEIGHING

# A Range of Veterinary Scales for Small and Large Animals





Cost-effective professional veterinary scales for general purpose weighing of pets and small animals.

## Small Animals



**MS-15:** The scale has a “HOLD” feature that will capture an unstable weight and keep it on the display for recording after the weight is taken off.

- **Capacity:** 44l b x 0.01 lb / 20 kg x 5 g
- **Construction:** ABS plastic with removable weighing tray
- **Inside Tray Dimension:** 23" x 11" / 584 mm x 292 mm
- **Power:** Four AA batteries (included)



**MS-20 Series:** The scale has a “HOLD” feature that will capture an unstable weight and keep it on the display for recording after the weight is taken off.

- **Capacity:** 44 lb x ½ oz / 20 kg x 10 g
- **Construction:** Stainless steel or plastic weighing tray. Plastic housed indicator with 16"-26" coiled cable.
- **Inside Tray Dimension:**  
MS20: 20 ½" x 10" x 2 3/8" / 521 mm x 254 mm x 60 mm  
MS20S: 20 ½" x 11 1/2" x 2 ¼" / 521 mm x 292 mm x 57 mm
- **Power:** 9V 100mAAC adapter or 9V alkaline battery (not included)

Ideal for accurately weighing large animals, these scales help to remove the guesswork where treatments or administration of feed needs to be closely monitored.

## Large Animals



**PS500-42:** This scale comes with a green rubber mat ideal for vet and other animal weighing applications.

- **Capacity:** 500 lb x 0.2 lb / 250 kg x 0.1 kg
- **Construction:** Low profile stainless steel platform with mild steel understructure and stainless steel/ABS indicator and a rubber slip mat
- **Platform Dimension:**  
42" x 22" x 2.7" / 1067 mm x 559 mm x 51 mm
- **Power:** Four AAA batteries or 9 VDC, 600mA adapter (both included)



**PS2000:** This scale is accurate, durable, reliable, portable and cost effective putting you at ease knowing you're in good hands.

- **Capacity:** 2000 lb x 1.0 lb / 1000 kg x 0.5 kg
- **Construction:** Rugged steel with tread plate surface and stainless steel/ABS indicator
- **Platform Dimension:**  
59" x 30" x 2.5"  
1498 mm x 762 mm x 64 mm
- **Power:** 9 VDC, 600mA adapter (included) or four AAA batteries (not included)

**PS1000:** This scale is accurate, durable, reliable, portable and cost effective putting you at ease knowing you're in good hands.



- **Capacity:** 1000 lb x 0.5 lb / 500 kg x 0.2 kg
- **Construction:** Rugged steel with tread plate surface and stainless steel/ABS indicator
- **Platform Dimension:**  
55.75" x 20.25" x 2.5" / 1416 mm x 514 mm x 64 mm
- **Power:** 9 VDC, 600mA adapter (included) or four AAA batteries (not included)

**PS3000-LCD:** This scale comes with a rubber mat that reduces noise and risk of animals/good slipping while weighing.



- **Capacity:** 3000 lb x 1.0 lb / 1500 kg x 0.5 kg
- **Construction:** Mild steel construction with safety tread plate surface
- **Platform Dimension:**  
78.75" x 39.37" x 3.93"  
2000 mm x 1000 mm x 100 mm
- **Power:** 9V 600mA universal power supply with USA/UK/EU adapters (included) or six AA batteries (not included)



## Laboratory Scales

**MBS Series:** With its fine feature set and simple functionality, this precision lab balance provides accurate weighing and counting solutions.



• **Capacity:**

300 g x 0.005 g, 600 g x 0.01 g,  
1200 g x 0.02 g, 3000 g x 0.05 g, 6000 g x 0.1 g

• **Construction:**

Robust plastic housing, stainless steel top plate

• **Platform Dimension:**

300 to 600 g: 4.65" x 0.4" / 120 mm x 9.5 mm  
1200 to 6000 g: 6.9" x 6.9" x 0.3" / 175 mm x 175 mm x 8 mm

• **Power:**

AA batteries (included) with auto-off

## Hanging Scales

**ElectroSamson:** This accurate hand-held scale featured with "Hold" mode and memory tare function is ideal for poultry weighing, seed and animal feed preparation.

• **Capacity:** 10 kg x 0.01 kg / 22 lb x 0.02 lb/¼ oz

25 kg x 0.02 kg / 55 lb x 0.05 lb/½ oz

45 kg x 0.05 kg / 99 lb x 0.1 lb/1 oz

• **Construction:** Tough ABS case with built-in hand grip for easy fitting

• **Product Dimension:**

11" x 5.5" x 3" / 279 mm x 140 mm x 76 mm

• **Power:** 9V alkaline battery (included)

