

 **Brecknell**

CS-200M

Chair Scale



The CS-200M is a mechanical chair scale designed for weighing those who have difficulty standing on their own. The die-cast easy-to-read beam runs along the back of the chair. Powder coated steel construction.

The chair is robust and oversized for comfort with armrests that lift away for easy access. The rear wheels lock for safety and the front wheels swivel for easy portability along with easy-to-clean plastic seat.

The CS-200M is the perfect solution for assisted living facilities, nursing homes and physician exam rooms.



Construction

Powder coated steel



Capacity

440 lb x 4 oz / 200 kg x 100 g



Display

Accurate full range weights that are dual marked in lb/oz and kg/g die-cast beam

VALUE IN
WEIGHING

CS-200M Chair Scale

Specifications

Capacity: 440 lb x 4 oz / 200 kg x 100 g

Readability: Accurate full range weights dual marked in lb/kg die-cast beam

Features: Locking rear wheels for safety. Swivel front wheels for transportation. Lift way armrests for easy access.

Seat Dimension (L x W x H):

16" x 14" x 13"

405 mm x 350 mm x 340 mm

Product Weight: 70 lb / 32 kg

Shipping Weight: 80 lb / 36 kg

Package Dimension (L x W x H):

30" x 24" x 25"

762 mm x 609 mm x 635 mm

Packaging Contents: Manual, footrest, two armrests, scale with wheel assembly, seat, back handle, weighing tower and miscellaneous hardware

Warranty: 3-year limited



 **Brecknell®**

brecknellscales.com

brecknellscales.co.uk

Part of Avery Weigh-Tronix

Brecknell is part of Avery Weigh-Tronix. Avery Weigh-Tronix is a trademark of Illinois Tool Works group of companies whose ultimate parent company is Illinois Tool Works Inc ("Illinois Tool Works:"). Copyright © 2023 Illinois Tool Works. All rights reserved. This publication is issued to provide outline information only and may not be regarded as a representation relating to the products or services concerned. This publication was correct at the time of going to print, however Avery Weigh-Tronix reserves the right to alter without notice the specification, design, price or conditions of supply of any product or service at any time.

CS-200M_L_500787.indd
V3 AWT35-500787