



The MS-15 scale is particularly suited for medical or veterinary applications.

With a capacity of 44 lb / 20 kg and readability of 0.01 lb / 5 g, the MS-15 will accurately weigh infants or small animals.

The scale has a "HOLD" feature that will capture an unstable weight and keep it on the display for recording after the weight is taken off.

The LCD display has a blue backlight with standard features of kg and lb and a hold feature.

Operates on four AA batteries (included). Easy-to-clean plastic weighing platter.



Accuracy

Within 0.2% of full load



Capacity

44 lb x 0.01 lb / 20 kg x 5 g



Power

Four AA batteries (included)



Display

1.0" / 25 mm LCD display with four active digits



Construction

ABS plastic with removable weighing tray



Operator Keys

On/Tare, Mode, Memory

MS-15 Medical/Veterinary Scale

Specifications

Unit: lb, kg

Weighing Range: 0.04 lb to 44 lb / 0.02 kg to 20 kg

Operating Temperature: 41° F to 95° F / 5° C to 35° C

Inside Tray Dimension (L x W):

23" x 11"

584 mm x 292 mm

Product Dimension (L x W):

23" x 13"

578 mm x 334 mm

Product Weight: 6.6 lb / 3 kg

Packaging Dimension (L x W x H):

25" x 15" x 6"

61 cm x 37 cm x 14 cm

Shipping Weight: 9 lb / 4 kg

Certifications: Full FCC compliance

Packaging Contents: Scale, batteries and manual

Warranty: 3-year limited



Brecknell®

brecknellscales.com

brecknellscales.co.uk



Brecknell is part of Avery Weigh-Tronix. Avery Weigh-Tronix is a trademark of Illinois Tool Works group of companies whose ultimate parent company is Illinois Tool Works Inc ("Illinois Tool Works:"). Copyright © 2023 Illinois Tool Works. All rights reserved. This publication is issued to provide outline information only and may not be regarded as a representation relating to the products or services concerned. This publication was correct at the time of going to print, however Avery Weigh-Tronix reserves the right to alter without notice the specification, design, price or conditions of supply of any product or service at any time.

Part of Avery Weigh-Tronix

MS-15_L_500801.indd
V4 AWT35-500801