


VALUE IN  
WEIGHING



A Range of Medical  
Scales for Home or  
Professional Use







Easily keep track of your infants weight at home or at the doctor's office with one of these pediatric scales. They quickly, safely and accurately weigh given you the peace of mind that your infant is at a healthy weight.

## Pediatric Scales

**MS-15:** This scale features a large LCD display and a removable top plate. The weight hold function freezes and stores the displayed weight.

- **Capacity:** 44 lb x 0.01 lb / 20 kg x 5 g
- **Construction:** ABS plastic with removable weighing tray
- **Tray Dimension:** 23" x 11" / 584 mm x 292 mm
- **Power:** Four AA batteries (included)
- **Unit:** lb, kg



MS-15 with  
removable top plate

**MS20 Series:** Comes with a plastic or a stainless steel weighing tray. The weight hold function freezes and stores the displayed weight.

- **Capacity:** 44 lb x ½ oz / 20 kg x 10 g
- **Construction:** Stainless steel or plastic weighing tray
- **Tray Dimension:** 21" x 10" / 533 mm x 254 mm
- **Power:** 9V 100mA AC adapter or 9V alkaline battery (not included)
- **Unit:** lb, kg



MS20

MS20-S

Choose from mechanical scale that can be placed anywhere or a digital scale that is easy-to-use and can transmit information.

## Portable Physician Scales

**HS-200M:** This mechanical scale provides you with the ability to take weight and height measurements.

- **Capacity:** 440 lb x 4 oz / 200 kg x 100 g
- **Construction:** Steel with durable white finish. Integrated wheel kit for easy portability.
- **Platform Dimension:**  
14.75" x 10.8" / 375 mm x 275 mm

**HS-300:** This scale provides accurate weight in lb/kg on easy-to-read LCD backlit display. Handrail and height rod are included.

- **Capacity:** 660 lb x 0.2 lb / 300 kg x 0.1 g
- **Construction:** Mild steel with durable powder coated white finish
- **Platform Dimension:**  
20.5" x 20.5" / 521 mm x 521 mm
- **Power:** Four AAA batteries or 9 VDC, 600mA power adapter (both included)

**MS140-300:** This digital handrail scale is the ideal solution for any application requiring accuracy and medium to heavy weight capacity.

- **Capacity:** 660 lb x 0.2 lb, 660 lb x 5 oz, 300 kg x 0.1 g
- **Construction:** Stainless steel weighing surface with mild steel foam covered handrails
- **Platform Dimension:**  
20.5" x 20.5" / 521 mm x 521 mm
- **Power:** Four AAA batteries or 9 VDC, 600mA power adapter (both included).



## Physician Scales

Bariatric and chair scales are ideal for any busy medical centers where the patient's weight needs to be regularly taken and monitored.

**MS1000-LCD:** This platform provides a safe and easy way to push patients in wheelchairs on and off the scale. Large easy-to-read LCD indicator and stainless steel handrail are included.

- **Capacity:** 1000 lb x 0.2 lb, 500 kg x 0.1 kg
- **Construction:** Powder coated rugged smooth steel deck structure and built-in ramps, fitted with two handles and wheels for relocating the scale
- **Platform Dimension:**  
45.7" x 41.3" / 1060 mm x 1050 mm
- **Power:** 9V 600mA power supply with US power adapter (included) or six AA batteries (not included)



**CS200-M:** This oversized mechanical chair scale is ideally suited for those who have difficulty standing on their own.

- **Capacity:** 440 lb x 4 oz, 200 kg x 0.1 kg
- **Construction:** Powder coated steel, easy-to-clean plastic seat
- **Seat Dimension:**  
16" x 14" x 13" / 405 mm x 350 mm x 340 mm



STAY  
CONNECTED  
WITH US!

